# Streamlined, Single-tube PCR Assay that Quantifies SMN1 and SMN2 Copy Numbers using Capillary Electrophoresis

<u>Walairat Laosinchai-Wolf</u>, Stela Filipovic-Sadic, Liangjing Chen, Shobha Gokul, Justin Janovsky, Vivian Le, Melissa Church, Michael Dodge, John Milligan, and Gary J Latham Asuragen, Inc., Austin, TX

# Summary

- Spinal Muscular Atrophy (SMA) is a genetic disease caused by deletion of or mutation in the SMN1 gene, where disease severity is modulated by the SMN2 copy number.
- We developed a prototype AmplideX® PCR/CE SMN1/2 kit<sup>†</sup> in a single-tube PCR assay, providing accurate quantification of SMN1 and SMN2 using capillary electrophoresis in less than 3 hours.
- Here we describe feasibility data demonstrating high specificity, broad input range, and reliable genotyping to more than 3 copies of both SMN1 and SMN2, including detection of rare hybrid genes (gene conversions).

# Introduction

Spinal Muscular Atrophy (SMA), an autosomal recessive neuromuscular disease caused by a loss of SMN1 gene function, is the primary genetic cause of infant death. The copy number of the highly similar SMN2 gene is an important predictor of the severity of SMA, as the SMN2 gene can produce some functional SMN protein that partially restores biological function. The antisense oligonucleotide nusinersen (marketed as SPINRAZA) promotes SMN2 alternative splicing to enhance the effectiveness of SMN2 functional replacement of SMN1. Early diagnosis of SMA, along with knowledge of SMN2 copy number, is critical for effective medical management. Herein, we report the performance of a prototype AmplideX PCR/CE SMN1/2 kit<sup>†</sup>, a single-tube PCR assay that quantifies SMN1 and SMN2 copy number including SMN1/2 hybrids (consistent with gene conversion) using capillary electrophoresis (CE).

# **Materials and Methods**

We developed a multiplexed PCR that simultaneously amplifies SMN1, SMN2, and an endogenous control (EC) in a single well. This PCR can also identify SMN1/2 gene conversion events. The PCR products were separated and quantified via Applied Biosystems™ 3500xl Genetic Analyzer with POP7 polymer with a 2.5kV, 20sec injection and 40min run time. The copy number of SMN1, SMN2, or SMN1/2 hybrid was calculated as the peak area ratio of target gene and EC normalized to a calibrator. If gene conversion was detected, SMN1 and SMN2 copy number was determined from the sum of SMN1 and an SMN1 hybrid and the sum of SMN2 and an SMN2 hybrid, respectively, indicating number of copies of exon 7 for each gene.



Figure 1. Assay Workflow. The workflow is streamlined with total assay time less than 3 hours. Total hands-on time is 45 minutes. CE instrument time is for a single injection, or 24 samples using an Applied Biosystems\* 3500xl Genetic Analyzer.

### Results



Figure 2. Assay Outputs Quantify SMN1 and SMN2 Copy Numbers Along with Hybrid Genes. A) Diagram illustrating gene conversion of exon 7 between SMN1 and SMN2. B) CE trace of sample with SMN1 Hybrid. C) CE trace of sample with SMN2 Hybrid. D) Examples of SMN1 and SMN2 copy numbers determined by the prototype assay using cell-line, blood and buccal DNA. SMN2>1 conversion was detected in blood #1 and #2 samples.




Figure 3. Specificity Study. Two DNA samples derived from Coriell cell lines with the absence of either SMN1 or SMN2 were tested with DNA inputs up to 100 ng. Expected SMN1 and SMN2 copy numbers for NA03813 are 0 and 3, respectively, and for NA17117 are 3 and 0. SMN1 and SMN2 copy number for both DNA samples determined by the assay agreed with expected copy numbers up to 100 ng DNA input. The assay is specific for SMN1 and SMN2 detection.




Figure 4. DNA Input Range. DNA input ranging from 10 to 100 ng was evaluated using two cell line samples (Coriell) and four blood samples. The expected SMN1 and SMN2 copy number was determined using independent methods.\(^1\) White bars designate bins for calling SMN1 and SMN2 copy number. All samples with DNA input from 10 to 100 ng produced expected copy numbers (ranging from 1 to  $\geq$  4).




Figure 5. Accuracy Study with 201 Samples. Samples derived from cell line, blood and buccal were tested and compared to an independent assay. White bars designate bins for calling SMN1 and SMN2 copy number. Expected copies were determined using an orthogonal method. The percent agreement between the two methods was 100% for both SMN1 and SMN2.

### Conclusions

- The prototype assay described here detects and accurately quantifies SMN1 and SMN2 copy numbers, including hybrid genes, within 3 hours including 45 minutes hands-on time.
- Copy numbers were accurately quantified using 10 to 100 ng DNA, and were reliably determined with DNA derived from >200 cell line, blood and buccal cell sources.
- The AmplideX PCR/CE SMN1/2 assay offers a rapid and robust single-tube PCR with reduced complexity compared to existing methods.

### Reference

1. Stabley DL, Harris AW, Holbrook J, et al. SMN1 and SMN2 copy numbers in cell lines derived from patients with spinal muscular atrophy as measured by array digital PCR

