

FOR REFERENCE PURPOSES

This manual is for Reference Purposes Only. DO NOT use this protocol to run your assays. Periodically, optimizations and revisions are made to the kit and protocol, so it is important to always use the protocol included with the kit.

NEXTflex™ 16S V4 Amplicon-Seq Kit 2.0 - 4

(Illumina Compatible)

Catalog #4203-01 (8 reactions)



This product is for research use only.

This manual is proprietary to Bioo Scientific Corp., and intended only for customer use in connection with the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose without the prior written consent of Bioo Scientific. Periodic optimizations and revisions are made to kit components and manuals. Follow the protocol included with the kit.

Bioo Scientific makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability for this product, or of the fitness of the product for any purpose. Bioo Scientific shall not be liable for any damages, including special or consequential damages, or expense arising directly or indirectly from the use of this product.

Bioo Scientific, NEXTflex, AIR, The NGS Experts, qRNA, and NanoQ are trademarks or registered trademarks of Bioo Scientific. All other brands and names contained herein are the property of their respective owners.

NEXTflex™ 16S V4 Amplicon-Seq Kit 2.0 - 4 - 4203-01

GENERAL INFORMATION	2
Product Overview.....	2
Contents, Storage and Shelf Life.....	2
Required Materials not Provided.....	3
Warnings and Precautions.....	3
NEXTflex™ 16S V4 AMPLICON-SEQ 2.0 PREPARATION PROTOCOL	4
NEXTflex™ 16S V4 Amplicon-Seq 2.0 Sample Preparation Flow Chart.....	4
Starting Material.....	4
Reagent Preparation.....	4
STEP A: PCR I Amplification.....	5
STEP B: PCR I Cleanup.....	6
STEP C: PCR II Amplification.....	7
STEP D: PCR II Cleanup.....	8
LIBRARY VALIDATION	9
APPENDIX A	10
Oligonucleotide Sequences.....	10
Reverse Primer Index Sequences and Reverse Complements.....	10
Low Level Multiplexing.....	10
RELATED PRODUCTS	11
NOTES	12

Product Overview

The NEXTflex™ 16S V4 Amplicon-Seq Kit 2.0 is designed to prepare multiplexed amplicon libraries that span the fourth hypervariable domain of microbial 16S ribosomal RNA (rRNA) genes. These libraries are compatible with paired-end sequencing on the Illumina® MiSeq platform.

There are two main steps involved in 16S V4 amplicon processing: an initial PCR amplification using customized PCR primers that target the V4 domain, and a subsequent PCR amplification that integrates relevant flow cell binding domains and unique 12 base pair sample indices. A limited number of cleanup steps ensures maximum recovery of amplicons for downstream sequencing. Custom sequencing read primers are not required with this version of the 16S V4 Amplicon-Seq Kit.

Contents, Storage and Shelf Life

The NEXTflex™ 16S V4 Amplicon-Seq Kit 2.0 contains enough material to prepare eight 16S V4 samples from genomic DNA for Illumina® compatible sequencing. The shelf life of all reagents is 12 months when stored properly. All components can be safely stored at -20°C.

Kit Contents	Amount
GREEN CAP	
NEXTflex™ Amplicon PCR Master Mix	240 µL
ORANGE CAP	
NEXTflex™ 16S V4 PCR I Primer Mix	16 µL
YELLOW CAP	
NEXTflex™ PCR II Barcoded Primer Mix 1 – 4	4 µL
WHITE CAP	
Resuspension Buffer	1 mL
Nuclease-free Water	1.5 mL

Required Materials not Provided

- 1 ng - 50 ng high-quality genomic DNA in up to 33 μ L nuclease-free water
- 96 well PCR Plate Non-skirted (Phenix Research, Cat # MPS-499) or similar
- Adhesive PCR Plate Seal (BioRad, Cat # MSB1001)
- Agencourt AMPure XP 5 mL (Beckman Coulter Genomics, Cat # A63880)
- Magnetic Stand -96 (Ambion, Cat # AM10027) or similar
- Thermocycler
- 2, 10, 20, 200 and 1000 μ L pipettes / multichannel pipettes
- Nuclease-free barrier pipette tips
- Vortex
- 80% Ethanol, freshly prepared (room temperature)

Warnings and Precautions


Bioo Scientific strongly recommends that you read the following warnings and precautions. Periodically, optimizations and revisions are made to the components and manual. Therefore, it is important to follow the protocol included with the kit. If you need further assistance, you may contact your local distributor or Bioo Scientific at nextgen@biooscientific.com.

- Do not use the kit past the expiration date.
- Ensure pipettes are properly calibrated as library preparations are highly sensitive to pipetting error.
- Try to maintain a laboratory temperature of 20°–25°C (68°–77°F).
- Genomic DNA sample quality may vary between preparations. It is the user's responsibility to utilize high quality Genomic DNA. Genomic DNA that is heavily nicked or damaged may cause library preparation failure. Absorbance measurements at 260 nm are commonly used to quantify DNA and 260 nm / 280 nm ratios of 1.8 - 2.0 usually indicate relatively pure DNA. Other quantification methods using fluorescent dyes may also be used. The user should be aware that contaminating RNA, nucleotides and single-stranded DNA may affect the amount of usable DNA in a sample preparation.
- It is required that NEXTflex™ 16S V4 PCR I & PCR II Primer Mixes are used during PCR amplification steps.

NEXTflex™ 16S V4 AMPLICON-SEQ 2.0 PREPARATION PROTOCOL

NEXTflex™ 16S V4 Amplicon-Seq 2.0 Sample Preparation Flow Chart

Figure 1: Sample flow chart with approximate times necessary for each step.



Starting Material

The NEXTflex™ 16S V4 Amplicon-Seq Kit 2.0 has been optimized and validated using 1 ng - 50 ng of high-quality bacterial genomic DNA.

Reagent Preparation

1. Briefly spin down each component to ensure material has not lodged in the cap or side of tube. Keep on ice and vortex each NEXTflex™ Mix just prior to use.
2. Allow Agencourt AMPure XP Beads to come to room temperature and vortex the beads until liquid appears homogenous before every use.

STEP A: PCR I Amplification

Materials

Bioo Scientific Supplied

GREEN CAP - NEXTflex™ Amplicon PCR Master Mix

ORANGE CAP - NEXTflex™ 16S V4 PCR I Primer Mix

WHITE CAP - Nuclease-Free Water

User Supplied

Thermocycler

96 Well PCR Plate

High-Quality Bacterial Genomic DNA, 1 ng - 50 ng

1. For each sample, combine the following reagents on ice in the PCR plate.

_ μL	High-Quality Bacterial Genomic DNA (1 ng - 50 ng in up to 33 μL)
_ μL	Nuclease-free Water
15 μL	NEXTflex™ Amplicon PCR Master Mix
2 μL	16S V4 PCR I Primer Mix
<hr/>	
50 μL	TOTAL

2. Mix reaction well by pipetting.
3. Apply adhesive PCR plate seal and place in thermocycler for the following PCR cycles:

4 min	95°	
30 sec	95°	
30 sec	56°	Repeat 10 cycles
90 sec	72°	
<hr/>		
4 min	72°	

STEP B: PCR I Cleanup

Materials

Bioo Scientific Supplied

WHITE CAP - Resuspension Buffer

User Supplied

Agencourt AMPure XP Magnetic Beads (room temperature)

80% Ethanol, freshly prepared (room temperature)

Magnetic Stand

1. Add 50 μL of AMPure XP Beads to each sample. Mix thoroughly by pipetting.
2. Incubate at room temperature for 5 minutes.
3. Place the 96 well PCR Plate on the magnetic stand at room temperature until the supernatant appears completely clear.
4. Remove and discard the supernatant. Do not disturb beads. Some liquid may remain in wells.
5. With plate on stand, add 200 μL of freshly prepared 80% ethanol to each magnetic bead pellet and incubate plate at room temperature for 30 seconds. Carefully, remove ethanol by pipette.
6. Repeat step 5, for a total of 2 ethanol washes. Ensure all ethanol has been removed.
7. Remove the plate from the magnetic stand and let dry at room temperature for 3 minutes.
8. Resuspend dried beads with 17 μL of Resuspension Buffer. Mix thoroughly by pipetting. Ensure beads are no longer attached to the side of the well.
9. Incubate resuspended beads at room temperature for 2 minutes.
10. Place plate on magnetic stand for 5 minutes or until the sample appears clear.
11. Transfer 16 μL of clear supernatant (purified PCR I product) to new well.

STEP C: PCR II Amplification

Materials

Bioo Scientific Supplied

GREEN CAP - NEXTflex™ Amplicon PCR Master Mix

YELLOW CAP - NEXTflex™ PCR II Barcoded Primer Mix

WHITE CAP - Nuclease-Free Water

User Supplied

Thermocycler

96 Well PCR Plate

Purified PCR I product (from STEP B)

1. For each sample, combine the following reagents on ice in the PCR plate.

16 µL Purified PCR I product (from STEP B)

17 µL Nuclease-free Water

15 µL NEXTflex™ Amplicon PCR Master Mix

2 µL NEXTflex™ PCR II Barcoded Primer Mix

50 µL TOTAL

2. Mix well by pipette.
3. Apply adhesive PCR plate seal and place in thermocycler for the following PCR cycles:

4 min 95°

30 sec 95°

30 sec 60° *Repeat cycles as recommended in table below*

30 sec 72°

4 min 72°

Input to PCR I (ng)	PCR II Cycles
1	22
5	20
10	18
25	14
50	12

STEP D: PCR II Cleanup

Materials

Bioo Scientific Supplied

WHITE CAP - Resuspension Buffer

User Supplied

Agencourt AMPure XP Magnetic Beads (room temperature)


80% Ethanol, freshly prepared (room temperature)

Magnetic Stand

1. Add 50 μL of AMPure XP Beads to each clear sample. Mix thoroughly by pipetting.
2. Incubate at room temperature for 5 minutes.
3. Place the 96 well PCR Plate on the magnetic stand at room temperature until the supernatant appears completely clear.
4. Remove and discard the supernatant. Do not disturb beads. Some liquid may remain in wells.
5. With plate on stand, add 200 μL of freshly prepared 80% ethanol to each magnetic bead pellet and incubate plate at room temperature for 30 seconds. Carefully, remove ethanol by pipette.
6. Repeat step 5, for a total of 2 ethanol washes. Ensure all ethanol has been removed.
7. Remove the plate from the magnetic stand and let dry at room temperature for 3 minutes.
8. Resuspend dried beads with 17 μL of Resuspension Buffer. Mix thoroughly by pipetting. Ensure beads are no longer attached to the side of the well.
9. Incubate resuspended beads at room temperature for 2 minutes.
10. Place plate on magnetic stand for 5 minutes until the sample appears clear.
11. Transfer 16 μL of clear supernatant to new well.
12. To ensure cluster generation, it is recommended that you quantify your library by gel or Agilent Bioanalyzer. To quantify by gel, load 2 μL of 6X Gel Loading Dye and 6-10 μL of PCR Product on a 2% low melt agarose gel + SYBR Gold.
13. Quantitate DNA library templates for optimal cluster density. Libraries generated with this kit are sequenced with standard Illumina sequencing primers.

LIBRARY VALIDATION

Figure 2. Gel validation of the NEXTflex™ 16S V4 Amplicon PCR product.



*Important note – Bacterial hypervariable regions vary in base composition and length. For community studies, expect bands that are ~450 bp.

Oligonucleotide Sequences

NEXTflex™ 16S V4 PCR I Primer Mix	
NEXTflex™	Sequence 5' → 3'
16S V4 Forward	GACGCTCTCCGATCTTATGGTAATTGTGTGCCAGCMGCCGCGGTAA
16S V4 Reverse	TGTGCTCTTCCGATCTAGTCAGTCAGCCGGACTACHVGGGTWTCTAAT

NEXTflex™ 16S V4 PCR II Barcoded Primer Mix	
NEXTflex™	Sequence 5' → 3'
PCR II Forward	AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT
PCR II Reverse	CAAGCAGAAGACGGCATACGAGATXXXXXXXXXXXXGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

XXXXXXXXXXXX denotes the index region of adapter. The index sequences and the respective reverse complement sequences contained in each adapter are listed below.

Reverse Primer Index Sequences and Reverse Complements

Barcoded Primer	Index Sequence (5' → 3')	Reverse Complement
1	GGCCGGCTAGAT	ATCTAGCCGGCC
2	AAGGAAGAGATA	TATCTCTTCCTT
3	GGACGGCATCTA	TAGATGCCGTCC
4	AAGGAAGGAGCG	CGCTCCTTCCTT

For an electronic list of the 16S V4 Primers visit [our webpage](#).

Low Level Multiplexing

Use the following reverse primer combinations for low level multiplexing in this kit:

Pool of 2: (Barcodes 1 & 2) OR (Barcodes 3 & 4)

Pool of 3: (Barcodes 1, 2, & 3), (Barcodes 1, 2, & 4), (Barcodes 1, 3, & 4), OR (Barcodes 2, 3, & 4)

RELATED PRODUCTS

Illumina Compatible RNA NGS Kits and Adapters

NEXTflex™ Rapid Directional RNA-Seq Kit

NEXTflex™ RNA-Seq Barcodes

NEXTflex-96™ RNA-Seq Barcodes

NEXTflex™ Rapid Directional qRNA-Seq™ Kit

NEXTflex™ Small RNA Sequencing Kit v2

NEXTflex™ Small RNA Barcode Primers

NEXTflex™ Poly(A) Beads

Illumina Compatible DNA NGS Kits and Adapters

NEXTflex™ 16S V4 Amplicon-Seq Kit

NEXTflex™ 16S V4 Amplicon-Seq Kit 2.0

NEXTflex™ 16S V1-V3 Amplicon-Seq Kit

NEXTflex™ 18S ITS Amplicon-Seq Kit

NEXTflex™ Rapid DNA-Seq Kit

NEXTflex™ Cell Free DNA-Seq Kit

NEXTflex™ DNA Barcodes

NEXTflex-96™ DNA Barcodes

NEXTflex-HT™ Barcodes

NEXTflex™ Dual-Indexed DNA Barcodes

NEXTflex™ Bisulfite-Seq Kit

NEXTflex™ Bisulfite-Seq Barcodes

NEXTflex™ Methyl-Seq 1 Kit

NEXTflex™ Msp 1

NEXTflex™ ChIP-Seq Kit

NEXTflex™ ChIP-Seq Barcodes

NEXTflex-96™ ChIP-Seq Barcodes

NEXTflex™ Pre-Capture Combo Kit

NEXTflex™ Rapid Pre-Capture Combo Kit

NEXTflex™ DNA Barcode Blockers

NEXTflex™ PCR-Free DNA Sequencing Kit

NEXTflex™ PCR-Free Barcodes

NOTES



WE WANT TO HEAR FROM YOU!

Your feedback is important to us. Tell us what you think of our kits by scanning the QR code or visiting our website at www.biooscientific.com/NGSfeedback.

We can't wait to hear from you!



BIOO SCIENTIFIC®

THE NGS EXPERTS™

innovación
tecnológica
para
laboratorio

Rafer

www.rafer.es

Barcelona

93 645 50 28
barcelona@rafer.es

Bilbao

94 499 85 80
bilbao@rafer.es

La Coruña

981 93 89 26
galicia@rafer.es

Madrid

91 365 15 70
madrid@rafer.es

Málaga

639 359 792
malaga@rafer.es

Sevilla

954 369 334
sevilla@rafer.es

Valencia

96 340 48 00
levant@rafer.es

Zaragoza

976 23 74 00
rafer@rafer.es

Lisboa

21 154 19 98
lisboa@rafer.es