


Wako

microRNA Research

The Pathway of Translational Silencing by microRNA



The microRNA matured in multiple steps in the cells is taken into a protein compound called an RNA-induced silencing complex (RISC), and it is known that it bonds with Argonaute (Ago) subfamily protein as its main component. Specific purification of RNA bound to protein of the Ago subfamily is considered to be an important technology for microRNA function determination, discrimination of new microRNA, and analysis of interaction between microRNA and mRNA. Wako has developed and launched high-quality anti-Ago2 monoclonal antibodies, microRNA Isolation kits, and high-efficiency microRNA Cloning kits, and the line-up includes research reagents with powerful support for microRNA research.

for Highly Efficient Adapter Ligation

microRNA Cloning Kit Wako



Wako Cat. No. 290-66501

for Effective 65°C Ligation of ssDNA, ssRNA and ssRNA-ssDNA

ssDNA Ligase, thermostable



Wako Cat. No. 298-65103; 292-65101

Immunoprecipitation method with anti-Ago2 mAb


microRNA Isolation Kit, Human / Mouse Ago2 Anti Human / Mouse Ago2, Monoclonal Antibody



microRNA Isolation Kit, Human Ago2
(Wako Cat. No. 292-66701)



microRNA Isolation Kit, Mouse Ago2
(Wako Cat. No. 292-67301)



Anti Mouse Ago2, Monoclonal Antibody
(Clone No. 2D4), immunostained image

NEW METHOD!!

microRNA "Specific" Purification

microRNA Isolation Kit, Human Ago2 (Wako Cat. No. 292-66701), Mouse Ago2 (Wako Cat. No. 292-67301)

microRNA Isolation Kit, Human/Mouse Ago2 can prepare high purified fractions of microRNA, which binds with Argonaute2 (Ago2) protein, based on immunoprecipitation method by using a high affinity monoclonal antibody against Ago2.


The purified microRNA fraction will contain very little contaminated degradation fragments of rRNA and tRNA. These kits will highly improve the microRNA cloning efficiency compared with that by using conventional microRNA purification method.

Features

- 1 Achieve Immunoprecipitation of internal Ago2 from human, mouse, rat, and hamster. Note)
- 2 microRNA bound to Ago2 protein can be isolated with high purity.
- 3 Few rRNA or tRNA decomposition products or impurities like smallRNA etc.
- 4 The isolated microRNA fraction can be used for cloning and micro arrays.

Note) Wako Cat. #292-66701 is for human cells and tissues and #292-67301 is for cells and tissues of mice and rats.


Outline of procedure



Kit contents (10 reactions)


- 1 Anti Human Ago2 Antibody Beads Solution 500 μ L \times 1 vial
- 2 Cell Lysis Solution 50 mL \times 1 vial
- 3 Elution Solution 500 μ L \times 1 vial
- 4 Ethachinmate 30 μ L \times 1 vial
- 5 3mol/L Sodium Acetate 400 μ L \times 1 vial

Simple procedure



Comparison with conventional method

	Conventional method	Wako's Latest Method
Total RNA extraction	Necessary	Unnecessary
small RNA (≤ 200 nt) purification		
Denaturing PAGE		
Gel extraction from denaturing PAGE		
Gel purification from denaturing PAGE		
Contamination ratio of rRNA fragments	High	Low
Contamination ratio of tRNA fragments	High	Low
Handling time	≤ 12 hr	≤ 4 hr



Specific Purification of microRNA fraction from several cell lines




Figure. 1 Purification of microRNA fraction by using *microRNA Isolation Kit, Human Ago2* (Wako catalog #292-66701).
The purified microRNA fraction from human cultured cell lines (HeLa, HepG2, HEK293) were specifically detected by Urea-PAGE. Cell number of each cell line is approximately 5×10^6 . The applied volume per lane is half of $10 \mu\text{L}$ of final solution prepared with an IP by this kit.




Figure. 2 Urea-PAGE pattern of purified RNA by using *microRNA Isolation Kit, Mouse Ago2* (Wako catalog #292-67301).
The purified microRNA fractions from cultured rodent cell lines (P388D1, CHO-K1, PC-12) were detected by silver stain. Cell number of each cell line is approximately 1×10^7 . The applied volume per lane is half of isolated sample by this kit.

Cloning of purified microRNA

High efficiency of microRNA cloning by using these kits followed by using *microRNA Cloning Kit Wako*.




Figure. 3 Cloning efficiency of microRNA from HeLa cell lysate
The microRNA fraction was prepared by *microRNA Isolation Kit, Human Ago2* (#292-66701) from 5×10^6 HeLa cells. The cDNA encoding microRNA was synthesized by *microRNA Cloning Kit Wako* (#290-66501) and was inserted into T-vector. The presence ratio of microRNA was more than 90%. The contents of cloned microRNA species were indicated on Table 1.




Figure. 4 Cloning efficiency of microRNA from P388D1 cell lysate
The microRNA fraction was prepared by *microRNA Isolation Kit, Mouse Ago2* (#292-67301) from 5×10^6 P388D1 cells. The cDNA encoding microRNA was synthesized by *microRNA Cloning Kit Wako* (#290-66501) and was inserted into T-vector. The presence ratio of microRNA was more than 90%. The contents of cloned microRNA are indicated on Table 2.

Purification of microRNA fraction from tissues




Figure. 5 microRNA fractions, which were isolated by this kit from human liver and testis lysate, were detected by silver stain after Urea-PAGE.

(50mg Tissues → 2mL Total tissue lysate → 1mL Tissue lysate)
(→ Immunoprecipitation → 10μL Total RNA solution → 5μL RNA solution → Urea-PAGE.)




Figure. 6 microRNA fractions, which were isolated by this kit from mouse lung, liver and brain lysate, were detected by silver stain after Urea-PAGE.

(50mg Tissues → 2mL Total tissue lysate → 1mL Tissue lysate)
(→ Immunoprecipitation → 10μL Total RNA solution → 5μL RNA solution → Urea-PAGE.)

GOOD METHOD!

Achieve Highly Efficient, Accurate & Easy-to-Use Adapter Ligation

microRNA Cloning Kit *Wako* (Wako Cat. No. 290-66501)

The **microRNA Cloning Kit *Wako*** can prepare the cDNA encoding microRNA. The cloning procedure will be completed within 1.5 days after preparation of microRNA fraction.

This kit is supported by shrimp alkaline phosphatase (SAP), single strand DNA Ligase, thermostable (Wako Cat. No. 298-65103; 292-65101) (selling separately), and original modified adaptors.

The cloning efficiency using this kit is improved higher than that of the conventional methods, which used bacterial alkaline phosphatase and T4 RNA ligase.

Features

- Highly efficient and accurate adapter ligation by thermostable ligase.
- Suitable for cloning of microRNA forming secondary structures.
- Simple operation achieves forming of cDNA coding microRNA in 1.5 days.

Kit contents (8 reactions)

① SAP	16 μL	⑩ 1mol/L Tris-HCl, pH7.5	160 μL
② 5×SAP Buffer	64 μL	⑪ Ethachinmate	24 μL
③ 40×Ligation Buffer	16 μL	⑫ 10mol/L Ammonium Acetate	960 μL
④ RNase Inhibitor	16 μL	⑬ 3' Adaptor (50pmol/μL)	8 μL
⑤ 10mmol/L MnCl ₂	16 μL	⑭ 5' Adaptor (50pmol/μL)	8 μL
⑥ Reverse Transcriptase	8 μL	⑮ RT Primer (50pmol/μL)	8 μL
⑦ 10×RT Buffer	16 μL	⑯ 5' PCR Primer (50pmol/μL)	16 μL
⑧ dNTP Mixture	112 μL	⑰ 3' PCR Primer (50pmol/μL)	16 μL
⑨ 0.5mol/L EDTA, pH8.0	16 μL	⑱ Control RNA (30ng/μL)	8 μL

Thermostable ligase and DNA polymerase are not included in this kit.


IMPORTANT

microRNA Cloning Kit *Wako* must be used together with **Single Strand DNA Ligase, Thermostable (Wako Cat. No. 298-65103; 292-65101)**, which is sold separately as a single product.

This ligase will be used for the high efficient adaptor ligation and reaction buffers of this kit are optimized with this ligase.

This ligase is ATP-dependent and used for ligation of ssDNA-ssDNA and ssRNA-ssRNA. Reaction of this ligase at 55-65°C, which is the optimal temperature of enzymatic reaction, will promote the ligation efficiency of single strand nucleic acids.

Outline of procedure



Comparison with conventional method

	Conventional method	microRNA Cloning kit <i>Wako</i>
Number of adaptor ligation steps	8	3
Phosphorylation of adaptors	Necessary	Unnecessary
Complete inactivation of Alkaline Phosphatase	Impossible	Possible
Apparatus for magnet beads	Necessary	Unnecessary
Rl	Necessary	Unnecessary
BioAnalyzer (Agilent)	Necessary	Unnecessary
Adaptor ligation time	≥ 8 hr	2.5 hr
Handling time	2.5 Days	1.5 Days
Reproducibility of adaptor ligation time	Low	High
Efficiency of microRNA cloning	< 10%	> 70%
Detection by EtBr	Impossible	Possible
Cloning efficiency of secondary structured microRNA	Low (by T4 RNA Ligase)	High (by Single Strand DNA Ligase, Thermostable)

ssDNA Ligase, thermostable (Wako Cat. No. 298-65103; 292-65101)



- 1 High thermal stability
- 2 Optimum temperature: 55~65°C
- 3 High ligation efficiency
- 4 Ligation of ssRNA, ssDNA, and ssRNA-ssDNA is applicable.

Source: *E. coli* expressed thermophilic phage TS2126 single strand DNA ligase
 Appearance: 10mmol/L Tris-HCl (pH 8.0), 50mmol/L KCl, 0.1mmol/L EDTA, 1mmol/L DTT and 50% Glycerol
 Activity: Shown on each label. (approx. 10units/μL)

Application with the microRNA Cloning Kit *Wako* is described in the package insert of the kit.

Cloning of microRNA from HeLa cells

Procedure

- 1 Preparation of total RNA from HeLa cells (1×10^7 cells) by ISOGEN (Nippon Gene #315-02504, 10mL).
- 2 Preparation of small RNA fraction, less than 200nt, from total RNA by microRNA Isolation Kit (Bio Chain Institute Inc. catalog #KS341025).
- 3 Separation of microRNA fraction by denaturing PAGE.
- 4 Collection of the gels of 20~23nt region after electrophoresis.
- 5 Cloning by using microRNA Cloning Kit *Wako* (Wako catalog #290-66501).
- 6 Construction of the plasmids harboring cDNA encoding microRNA and transformation of *E. coli*.
- 7 Random selection of the 96 transformed *E. coli* from selection LB agar medium.
- 8 Determination and verification of the cDNA sequences by using Sanger miRBase.





Table 3. Contents of cloned microRNA species

microRNA species	The number of clone
hsa-miR-23a	38
hsa-miR-92a	16
hsa-miR-22	5
hsa-miR-25	5
hsa-miR-23b	3
hsa-miR-19b	1
hsa-miR-21	1
hsa-miR-210	1
Total	70


Figure. 7 Cloning efficiency of microRNA from HeLa cell lysate

The cloning efficiency of microRNA was more than 70%. Others indicate that isolated cDNA sequences were not matched miRBase. Unknowns indicate that isolated cDNA sequences were not matched human genome sequence. The contents of cloned microRNA species are indicated on Table 3.

HIGH QUALITY!!

Anti Ago2, Monoclonal Antibody

Anti Human Ago2, Monoclonal Antibody (Clone No.4G8) (Wako Cat. No. 011-22033; 015-22031)




Lane M : Molecular weight marker Lane 3 : HEK293 (Human)
Lane 1 : HeLa (Human) Lane 4 : THP-1 (Human)
Lane 2 : HepG2 (Human) Lane 5 : P388D1 (Mouse)

Application	Working Dilution
Western Blot	1 : 100 - 1 : 200
Immunoprecipitation	1 : 50
Immunocytochemical staining	1 : 20 - 1 : 50

Species	Human	Mouse	Hamster	Rat
Cell	HeLa	P388D1	CHO	SCC-131
Western Blot	○	×	×	×
Immunoprecipitation	○	×	×	×
Immunocytochemistry	○	×	×	×
microRNA Purification	○	×	×	×

Figure. 8 Immunoprecipitation of hAgo2 protein from human cultured cell lines (HeLa, HepG2, HEK293, THP-1) and mouse cultured cell line (P388D1) by using 20μL 10% Protein G slurry immobilized with 10μg of this antibody. The bands of hAgo2 protein were detected in approximately 100kDa by using silver staining and western blot. Cell number was 5×10⁶ cells.

Anti Mouse Ago2, Monoclonal Antibody (Clone No.2D4) (Wako Cat. No. 014-22023; 018-22021)



Lane 1 : NIH-3T3 (Mouse)
Lane 2 : SCC-131 (Rat)
Lane 3 : CHO (Hamster)

Application	Working Dilution
Western Blot	1 : 200 - 1 : 1,000
Immunoprecipitation	10μg / IP
Immunocytochemical staining	1 : 100 - 1 : 500

Species	Mouse	Hamster	Rat	Human
Cell	P388D1 NIH-3T3	CHO	SCC-131	NCI-H460
Western Blot	○	○	○	×
Immunoprecipitation	○	○	○	×
Immunocytochemistry	○ (NIH-3T3)	NT	NT	×
microRNA Purification	○ (P388D1)	○	○	×

NT : Not Tested

Figure. 9 Immunoprecipitation of Ago2 protein from NIH-3T3(Mouse), SCC-131(Rat) and CHO(Hamster) cell line by using 20μL of 10% Protein G slurry immobilized with 5μg of this antibody (2D4). The bands of endogenous Ago2 protein were detected in approximately 100kDa by using silver staining and western blot. The 1/1,000 diluted antibody was used as the 1st antibody for western blot. Cell number was 5×10⁶ cells.

Product List

Description	Wako Cat. No. (Package Size)
NEW!! microRNA Isolation Kit, Mouse Ago2	292-67301 (10 reactions)
microRNA Isolation Kit, Human Ago2	292-66701 (10 reactions)
NEW!! Anti Mouse Ago2, Monoclonal Antibody	014-22023 (50 μL)
NEW!! Anti Mouse Ago2, Monoclonal Antibody	018-22021 (100 μL)
NEW!! Anti Human Ago2, Monoclonal Antibody	011-22033 (50 μL)
NEW!! Anti Human Ago2, Monoclonal Antibody	015-22031 (100 μL)
microRNA Cloning Kit Wako	290-66501 (8 reactions)
Single Strand DNA Ligase, thermostable, recombinant, Solution	298-65103 (200 units)
	292-65101 (500 units)

- microRNA Isolation Kit, Human/Mouse Ago2 is patent pending (11, 30, 2007).
- microRNA Cloning Kit **Wako** is patent pending (1, 10, 2007).
- All products shown on the leaflet are for research use only.

■ Please visit our online catalog to search for other products from Wako ; <http://www.e-reagent.com>

Wako Pure Chemical Industries, Ltd.
<http://www.wako-chem.com>
 1-2, Doshomachi 3-Chome
 Chuo-Ku, Osaka 540-8605, Japan
 Tel: 81-6-6203-3741
 Fax: 81-6-6203-1999

Wako Chemicals USA, Inc.
<http://www.wakousa.com>
Head Office(Richmond, VA):
 Toll-Free(U.S. only):1-877-714-1920
 Tel: 1-804-714-1920/Fax: 1-804-271-7791
Los Angeles Sales Office(Irvine, CA):
 Tel: 1-949-679-1700/ Fax: 1-949-679-1701
Boston Sales Office(Cambridge, MA):
 Tel: 1-617-354-6772 Fax: 1-617-354-6774

Wako Chemicals GmbH
<http://www.wako-chemicals.de>
European Office:
 Fuggerstraße 12, D-41468
 Neuss, Germany
 Tel: 49-2131-311-0
 Fax: 49-2131-311100