

Gene Conversions and Hybrid Peak Detection in AmplideX® PCR/CE SMN1/2 Kit*

Overview of the Kit

The AmplideX® PCR/CE SMN1/2 Kit is an *in vitro* nucleic acid amplification kit for the determination of SMN1 and SMN2 exon 7 copy number¹. It is a single-tube PCR assay that amplifies distinctive SMN1 and SMN2 gene regions and an endogenous control (EC) gene from purified genomic DNA.

In this method, fluorescently-labeled SMN1- and SMN2-specific amplicons are resolved by capillary electrophoresis (CE) and referenced to co-amplified EC gene products to determine copy numbers. Discrete amplicon peaks are generated from the EC, SMN1, SMN2, and/or SMN1 or SMN2 hybrid, or chimeric, genes. Example electropherograms are shown in Figure 1.

SMN1/2 Gene Conversions and Hybrid Peaks

Survival motor neuron genes (telomeric SMN1 and centromeric SMN2) are located within a complex region of chromosome 5q13, which is prone to duplications and deletions. These two genes have equivalent promoters and differ only by five nucleotides within their 3' ends (in intron 6, exon 7, intron 7, and noncoding exon 8)². Despite their high homology, only the SMN1 gene is necessary for survival of motor neurons. A single nucleotide difference (c. 840C>T) between SMN1 and SMN2 in exon 7 disrupts an exonic splicing enhancer and reduces full-length SMN2 transcripts³. Due to the functional impact of this critical difference when SMN1 is deleted, methods for dosage analysis of the SMN1 and SMN2 genes target copy number quantification of exon 7.




Figure 1. A) Commonly observed sample genotype with 2 copies of SMN1 and 2 copies of SMN2 gene (no hybrid peaks). B) Sample with SMN1 hybrid peak (SMN1 sequence in exon 7 with SMN2 in intron 7). C) Sample with SMN2 hybrid peak (SMN2 sequence in exon 7 and SMN1 in intron 7).

In addition to SMN1 and SMN2 gene dosage, the AmplideX PCR/CE SMN1/2 Kit is also designed to detect and differentiate SMN1-to-SMN2 and SMN2-to-SMN1 gene conversions^{4,5,6,7,12} within exon 7 and intron 7 regions. When present, these chimeras are detected as unique peaks in the CE trace, referred to as hybrid peaks.

An SMN1 hybrid peak indicates a chimeric gene with SMN1 sequence in exon 7 and SMN2 sequence in intron 7. Similarly, an SMN2 hybrid peak stems from SMN1 gene conversion resulting in a SMN2 sequence in exon 7 and an SMN1 sequence in intron 7 (Figure 2). Importantly, only exon 7 status is used in the final copy number calculation. The presence of exon 7/intron 7 sequence mismatches due to gene conversion (as indicated by hybrid peaks) is currently provided for informational use only.

AmplideX®
PCR/CE SMN1/2 Kit

The AmplideX PCR/CE *SMN1/2* Macro enables streamlined analysis of *SMN1* and *SMN2* copy numbers, and reports genotypes as 0, 1, 2, 3, or ≥ 4 copies of exon 7 for both *SMN1* and *SMN2*. The exon 7 copy number represents the sum of gene-specific and hybrid gene integer copy numbers for both *SMN1* and *SMN2* (Figure 1). In addition to *SMN1* and *SMN2* genotypes, the macro output indicates a gene conversion event within the exon 7/intron 7 region in a separate hybrid peak column (labeled as “*SMN1*,” “*SMN2*,” or “both” when gene conversion is present)⁸.

Frequency of *SMN1/2* Gene Conversion Events

Conventional qPCR-based *SMN1/2* assays query c.840C>T in exon 7 and cannot flag *SMN* hybrid genes. These chimeras may be detected if a method is used to interrogate multiple loci that distinguish the two genes, as is the case for the AmplideX PCR/CE *SMN1/2* Kit. We note, however, that the junctions for *SMN1/2* chimeras are heterogeneous⁷, and thus the identification of *SMN* hybrids across different assay designs that probe multiple gene-specific sequences may not agree with one another, even though each may produce an analytically accurate result.

Several studies have reported that *SMN1/2* gene conversions are present in approximately 1% of the population^{9,10,11}. A similar rate of gene conversion was detected using the AmplideX PCR/CE *SMN1/2* Kit; a screen of 1426 human non-SMA samples uncovered 14 samples with *SMN* gene conversions (0.98%). Ten of the conversions (0.70%) had *SMN1* exon 7 along with *SMN2* intron 7 (*SMN1* hybrid peak) and four samples (0.28%) had *SMN2* exon 7 with *SMN1* intron 7 (*SMN2* hybrid peak). Dosage analysis of hybrid peaks from this sample set showed that the majority of samples had only 1 copy of a hybrid peak but cases of 2 copies have also been detected (two out of the fourteen samples with gene conversions had 2 copies of a hybrid peak). The aggregate exon 7 *SMN1* and *SMN2* copy number results for all gene conversion samples were confirmed with a qPCR-based orthogonal method.

References

1. AmplideX® PCR/CE *SMN1/2* Kit Protocol Guide (Version 00002356v1)
2. Identification and characterization of a spinal muscular atrophy-determining gene. Lefebvre et al. Cell. 1995 Jan 13;80(1):155-65. PMID: 7813012
3. Technical standards and guidelines for spinal muscular atrophy testing. Prior TW, Nagan N, Sugarman EA, Batish SD, Braastad C. Genet Med. 2011 Jul;13(7):686-94. doi: 10.1097/GIM.0b013e318220d523. PMID: 21673580
4. Cure SMA website: www.cureSMA.org
5. Telomeric Region of the Spinal Muscular Atrophy Locus Is Susceptible to Structural Variations. Noguchi et al. Pediatr Neurol. 2016 May; 58:83-9. doi: 10.1016/j.pediatrneurol.2016.01.019. Epub 2016 Jan 30. PMID: 27268759
6. Inverse correlation between *SMN1* and *SMN2* copy numbers: evidence for gene conversion from *SMN2* to *SMN1*. Ogino S, Gao S, Leonard DG, Paessler M, Wilson RB. Eur J Hum Genet. 2003 Sep;11(9):723. PMID: 12939662
7. A new method for *SMN1* and hybrid *SMN* gene analysis in spinal muscular atrophy using long-range PCR followed by sequencing. Kubo Y, Nishio H, Saito K. J Hum Genet. 2015 May; 60(5):233-9. doi: 10.1038/jhg.2015.16. Epub 2015 Feb 26. PMID: 25716911
8. AmplideX® PCR/CE *SMN1/2* Macro User Guide (Version ASGN-SMN12-00)
9. Identification of bidirectional gene conversion between *SMN1* and *SMN2* by simultaneous analysis of *SMN* dosage and hybrid genes in a Chinese population. Chen TH, Tzeng C, Wang CC, Wu SM, Chang JG, Yang SN, Hung CH, Jong YJ. J Neurol Sci. 2011 Sep 15;308(1-2):83-7. doi:10.1016/j.jns.2011.06.002. Epub 2011 Jun 25. PMID: 21705024
10. Universal fluorescent multiplex PCR and capillary electrophoresis for evaluation of gene conversion between *SMN1* and *SMN2* in spinal muscular atrophy. Wang CC, Jong YJ, Chang JG, Chen YL, Wu SM. Anal Bioanal Chem. 2010 Jul;397(6):2375-83. doi: 10.1007/s00216-010-3761-1. Epub 2010 Jun 19. PMID: 20563565
11. SMA prenatal diagnosis: a modified protocol to help differentiation between deletions and gene conversion. Kekou K, Sofocleous C, Konstantinidis G, Fryssira H, Mavrou A, Kitsiou S, Kanavakis E. Mol Cell Probes. 2015 Feb;29(1):71-3. doi: 10.1016/j.mcp.2014.10.001. Epub 2014 Oct 13. PMID: 25308401
12. Complete sequencing of the *SMN2* gene in SMA patients detects *SMN* gene deletion junctions and variants in *SMN2* that modify the SMA phenotype. Ruhno C, McGovern VL, Avenarius MR, Snyder PJ, Prior TW, Nery FC, Muhtaseb A, Roggenbuck JS, Kissel JT, Sansone VA, Siranosian JJ, Johnstone AJ, Nwe PH, Zhang RZ, Swoboda KJ, Burghes AHM. Hum Genet. 2019 Feb 20. doi: 10.1007/s00439-019-01983-0. [Epub ahead of print] PMID: 30788592




Figure 2. Schematic of all *SMN* peak sizes in base pairs, and their corresponding exon 7/intron 7 status.