


1

Prepare Sample

- Pretreat the sample and extract the nucleic acid.
- Program the PCR protocol into the thermal cycler with a heated lid (105°C).
- Pre-heat the thermal cycler to 42°C.
- Setup the NxTAG® Gastrointestinal Pathogen Panel (GPP) reaction plate.
 - If frozen, thaw the extracted nucleic acid samples.
 - Vortex the samples and spin down.
 - Aspirate 35 µL of the sample or control to each PCR vessel.
 - Reseal the plate.



2

Run Thermal Protocol

- Place the foil-sealed plate in the pre-heated thermal cycler and run the protocol.
- Remove the plate from the thermal cycler.

3

Setup System Software

- Perform the following, if necessary:
 - Install the SYNCT™ Software.
 - Install the NxTAG module in SYNCT Software.
 - Import the NxTAG GPP assay file into SYNCT Software.
 - Import the NxTAG GPP data acquisition protocol into xPONENT® Software.
- Define the assay controls:
 - Negative amplification control
 - Negative control
 - Positive external control
- If allowed, define the Test Panel.

4

Acquire Samples

- In xPONENT, perform the **Enhanced Startup Routine** with probe sonication and adjust the sample probe height at least once a week.
- Heat the MAGPIX® plate heater to 37°C.
- Add the post-batch routine reagents.
- Place the plate on the heater block.
- Create a batch.
- Run the batch.
- Complete the run.
 - Turn off the plate heater and remove the plate from the heater block.
 - Discard the test strip(s) in a sealed biohazard bag.
 - Decontaminate the plate, if you are re-using it.
 - Run the **System Shutdown Routine**.

5

Analyze Results in SYNCT™ Software

- Create a run and import data into SYNCT Software.
- Edit and review the Orders.
 - For samples, choose a Test Panel OR customize by selecting or masking targets.
 - For controls, apply a pre-defined control.
- Process the Run by selecting the Sample ID (Run) and clicking **Process Run**.
- View the results by navigating to **Results > Results**.
- Create and print a report.

For Research Use Only. Not for use in diagnostic procedures.
NxTAG® GPP Quick Reference Guide (RUO)

MLD-053-MAN-001 Rev A 03/2021

© 2021 Luminex Corporation. All rights reserved.