Wako

For Genetic Research
MagCapture™ Exosome Isolation Kit PS

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Exosomes and other extracellular vesicles (EVs) are small membrane vesicles containing protein, mRNA, microRNA, DNA, and lipids, which are secreted by various cells and are stable in body fluids including blood, saliva, urine, cerebrospinal fluid (CSF) and breast milk. These EVs have been recognized as messengers of cell-to-cell communication, and biomarkers for various diseases. Conventionally, ultracentrifugation, affinity purification using antibody to surface antigen, and precipitation with polymer reagents are used for isolation of exosomes. However, these methods are not fully satisfying in all performances criteria such as recovery efficiency, purity and operability.

MagCapture™ Exosome Isolation Kit PS adopts a novel affinity purification method using magnetic beads and phosphatidylserine (PS)-binding protein (PS affinity method). This kit can easily isolate high purity exosomes and other EVs from cell culture medium and body fluids at high yield by a normal microcentrifuge. If higher purity exosomes are needed, please use the $10,000 \times g$ supernatant as sample. This kit enables the isolation of exosomes and other EVs as intact forms because the captured EVs are eluted from magnetic beads with metal-chelating reagent at neutral pH. The isolated intact exosomes and other EVs can be used for various applications including electron microscopic analysis, nanoparticle tracking analysis, administration of EVs and analysis of molecular constituents such as proteins, lipids, or nucleic acids.

【Features】

- High purity exosomes and other EVs are isolated by a novel affinity purification method (PS affinity method).
- The purity and yield of exosomes are higher than ultracentrifugation.
- This kit enables the purification of exosomes and other EVs as intact forms, and can be used for various applications.
- Exosomes can be isolated with magnetic beads.
- No ultracentrifugation is required.

【Kit contents】 (10 purifications)

This kit includes 7 components.

(1) Streptavidin Magnetic Beads	$600 \mu\text{L} \times 1$ tube
(2) Biotin-labeled Exosome Capture	$100 \mu\text{L} \times 1$ tube
(3) Exosome Capture Immobilizing Buffer	$35\text{mL} \times 1$ bottle
(4) Exosome Binding Enhancer ($\times 500$)	$500 \mu\text{L} \times 1$ tube
(5) Washing Buffer	$75\text{mL} \times 2$ bottles
(6) Exosome Elution Buffer	$5\text{mL} \times 1$ bottle
(7) Reaction Tubes	22 tubes

【Storage】

Store at 2-10 °C

【Additional required materials】

Reagents :

- 1) TBS (As necessary)

Equipment :

- 1) Microcentrifuge (Max $> 10,000 \times g$)
- 2) Vortex mixer
- 3) Tabletop centrifuge
- 4) Magnetic stand
- 5) Rotator or tube mixer
- 6) Micro pipette
- 7) Pipette tip
- 8) Microcentrifuge tube (1.5mL)
- 9) Centrifuge tube (15mL) (As necessary)
- 10) Centrifuge tube (50mL) (As necessary)
- 11) Ultrafiltration unit (Vivaspin6, M. W. cut off : 100K, Code No. VS0641 or Vivaspin20, M. W. cut off : 100K, Code No. VS2041) (As necessary)

【Precaution for use】

1. Equipment

Please use sterile or DNase and RNase free microcentrifuge tubes and pipette tips. We recommend the use of gloves and a mask to avoid contamination of DNase and RNase.


2. Reagents

If you analyze the nucleic acid (RNA, DNA) after isolation of exosomes and other EVs, please handle reagents carefully and keep them as sterile as possible. Additionally, if you analyze the mass spectrometry of proteins after isolation of exosomes and other EVs, please be careful of contamination from other proteins.

3. Handling of biohazardous waste

Please treat human serum, human plasma, and human tissue samples as infectious samples. Biohazardous wastes, which include waste liquid and equipment such as microcentrifuge tubes, pipette tips, and gloves, must be disposed of according to the guidelines of the institution.

【Outline of procedure】

【Procedure】

1. Preparation of sample

This is the section to prepare samples. If exosomes and other large EVs are needed, prepare **$1,200 \times g$ supernatant** as a sample. Additionally, if higher purity of exosomes are needed, prepare **$10,000 \times g$ supernatant** as a sample. This protocol for sample preparation is set for cell culture medium and serum. If other body fluids are used, please examine the pretreatment protocol by referencing the protocol for serum.

In the case of cell culture medium

- 1) Culture the cells under the appropriate condition^{※1}.
- 2) Harvest cell culture medium.
- 3) Centrifuge the cell culture medium for 5 minutes at $300 \times g$ and 4°C to remove cells and debris.
- 4) Transfer the supernatant from step 3) into a new tube.
- 5) Centrifuge 4) at $1,200 \times g$ for 20 minutes at 4°C to remove the cell debris.
- 6) Transfer the supernatant from step 5) into a new tube. (**$1,200 \times g$ supernatant**)
- 7) Centrifuge 6) at $10,000 \times g$ for 30 minutes at 4°C to remove the large EVs.
- 8) Transfer the supernatant from step 7) into a new tube. (**$10,000 \times g$ supernatant**)

[Optional]

If the volume of **$1,200 \times g$ supernatant** or **$10,000 \times g$ supernatant** is over 1mL, concentrating the sample by using ultrafiltration unit (Vivaspin6, M. W. cut off: 100K, Code No. VS0641 or Vivaspin20, M. W. cut off: 100K, Code No. VS2041) is recommended. Please concentrate until the sample is less than or equal to 1mL according to the manufacturer's instruction manual. The maximum sample volume of ultrafiltration treatment is 50mL (50-fold concentration). Exosomes and other EVs are efficiently recovered from concentrated sample.

※1 Please refer to appropriate culture conditions according to cell lines. To ensure that isolated exosomes originate from your cells, please use exosome-depleted FBS. Exosome-depleted FBS is prepared by ultracentrifugation (Example: at $110,000 \times g$ for 18 hours) or is commercially available.

In the case of serum


- 1) Centrifuge the serum^{※2} at $1,200 \times g$ for 20 minutes and 4°C to remove the cell debris.
- 2) Transfer the supernatant from step 1) into a new tube. (**$1,200 \times g$ supernatant**^{※3})
- 3) Centrifuge 2) at $10,000 \times g$ for 30 minutes at 4°C to remove the large EVs.
- 4) Transfer the supernatant from step 3) into a new tube. (**$10,000 \times g$ supernatant**^{※3})

※2 This kit is unsuitable for isolation from plasma treated with chelating agents (EDTA, citric acid).


※3 In the case of serum, $1,200 \times g$ supernatant and $10,000 \times g$ supernatant are unsuitable for concentration by using ultrafiltration unit.

【Flowchart of Sample Preparation (Section 1)】

Cell Culture Medium

Serum

2. Immobilization of Exosome Capture

This is the section to immobilize Biotin-labeled Exosome Capture reagent to Streptavidin Magnetic Beads. Please use the 1.5mL Reaction Tube (included) in the downstream experiment which requires magnetic beads operation^{※4}.


- 1) Transfer $60 \mu\text{L}$ ^{※5} of Streptavidin Magnetic Beads into a new 1.5mL Reaction Tube.
- 2) Place the 1.5mL Reaction Tube on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube.
- 3) Add $500 \mu\text{L}$ of Exosome Capture Immobilizing Buffer into the 1.5mL Reaction Tube in step 2) and suspend it by vortexing.
- 4) After spin-down of the 1.5mL Reaction Tube from step 3), place it on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube.
- 5) Add $500 \mu\text{L}$ of Exosome Capture Immobilizing Buffer and $10 \mu\text{L}$ of Biotin-labeled Exosome Capture into the 1.5mL Reaction Tube from step 4) and mix by vortexing.

- 6) Mix for 10 minutes at 2-10 °C with rotator or tube mixer.
- 7) After spin-down of the 1.5mL Reaction Tube from step 6), place it on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube.
- 8) Add 500 μ L of Exosome Capture Immobilizing Buffer into the 1.5mL Reaction Tube from step 7) and mix by vortexing.
- 9) After spin-down of the 1.5mL Reaction Tube from step 8), place it on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube.
- 10) Repeat steps 8)-9) twice. ··· **Exosome Capture-immobilized beads**

※4 Please only use the 1.5mL Reaction tube included in this kit because only these tubes can generate the high recovery rate of extracellular vesicles by preventing non-specific adsorption of magnetic beads to the tube wall.

※5 Please adequately suspend the magnetic beads by vortexing. Make sure homogeneous suspension of the magnetic beads is achieved before use.

[Flowchart of immobilization of Exosome Capture (Section 2)]

3. Affinity Reaction

This is the affinity reaction section of the Exosome Capture-immobilized beads with exosomes and other EVs in sample. This protocol is based on the use of 1.5mL Reaction Tube (max volume: 1mL). If exosomes and large EVs are needed, please use **1,200 × g supernatant** as sample. Additionally, if higher purity of exosomes are needed, please use **10,000 × g supernatant** as sample.

- 1) Transfer sample (max volume: 1mL) to a new 1.5mL microcentrifuge tube (not included in this kit) and add 1:500 volume of Exosome Binding Enhancer ($\times 500$) to sample. Mix by vortexing.
- 2) After spin-down of the 1.5mL microcentrifuge tube from step 1), transfer sample into the 1.5mL Reaction Tube containing **Exosome Capture-immobilized beads**. Mix by vortexing.
- 3) Mix for 3 hours~ at 2-10 °C with rotator or tube mixer.
- 4) After spin-down of the 1.5mL Reaction Tube from step 3), place it on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant^{※6} after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube. ··· **EVs-binding beads**

[Optional]

When the sample volume is over 1mL, please refer to the following protocol to use 15mL centrifuge tube (max volume: 10mL). However, the increase of sample volume may decrease recovery rate. Therefore, we recommend the use of the Exosome Capture-immobilized beads repeatedly in order to recover extracellular vesicles from used sample (refer to ※14). Please repeat sections 3-5 (up to 5 times). Additionally, in the case of cell culture medium, we recommend concentration less than or equal to 1mL by using ultrafiltration unit according to above-described protocol.

- 1) Transfer sample (max volume: 10mL) to a new 15mL centrifuge tube^{※7} (not included in this kit) and add 1:500 volume of Exosome Binding Enhancer ($\times 500$) to sample. Mix by vortexing.
- 2) Transfer 1mL of sample of 1) into the 1.5mL Reaction Tube^{※8} containing **Exosome Capture-immobilized beads**. After suspending the beads, transfer all of the beads in suspension into the 15mL centrifuge tube from 1). Mix by vortexing.
- 3) Mix for 3 hours~ at 2-10 °C with rotator or tube mixer.
- 4) Centrifuge 3) at 1,000 \times g for 3 minutes at 4 °C.
- 5) Remove the supernatant^{※6} carefully and leave 1ml of the supernatant and bead pellets.
- 6) After suspending sample from step 5), transfer all of the beads in suspension into the 1.5mL Reaction Tube from 2).
- 7) After spin-down of the 1.5mL Reaction Tube from step 6), place it on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant^{※6} after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube. ··· **EVs-binding beads**

※6 If further recovery of remaining extracellular vesicles is needed, please store the used supernatant.

※7 Please use a 15mL centrifuge tube in which the non-specific adsorption of magnetic beads to the tube wall is low. (Example: Code No. 430791, Corning)

※8 Please store this 1.5mL Reaction Tube after transferring the beads into the 15mL centrifuge tube, because it is used again at the optional step 6) of section 3.

4. Washing of EVs-binding beads

This is the washing section of **EVs-binding beads**. Please use Washing Buffer after adding Exosome Binding Enhancer ($\times 500$) certainly^{※9}.

- 1) Add 1:500 volume of Exosome Binding Enhancer ($\times 500$) to the required Washing Buffer^{※10} and mix by vortexing^{※11}.
 - 2) Add 1mL of Washing Buffer (+ Exosome Binding Enhancer ($\times 500$)) into the 1.5mL Reaction Tube containing **EVs-binding beads** and mix by vortexing.
 - 3) After spin-down of the 1.5mL Reaction Tube from step 1), place it on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube.
 - 4) Repeat steps 2)-3) twice.
 - 5) After spin-down of the 1.5mL Reaction Tube from step 4), place it on the magnetic stand for 1 minute to collect the beads. Remove and discard the supernatant completely after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube^{※12} . . .

Washed EVs-binding beads

- ※9 If the addition of Exosome Binding Enhancer ($\times 500$) to Washing buffer is omitted, captured exosome and EVs will be eluted from the magnetic beads by Washing Buffer during washing section.
 - ※10 The required volume of Washing Buffer is 3mL per reaction.
 - ※11 Please prepare Washing Buffer (+ Exosome Binding Enhancer ($\times 500$)) just before use, because constituents may be precipitated after the addition of Exosome Binding Enhancer ($\times 500$) overtime.
 - ※12 Please remove Washing Buffer (+ Exosome Binding Enhancer ($\times 500$)) in the tube completely, because remaining buffer causes reduction of elution efficiency of captured exosomes and other EVs.

5. Elution of extracellular vesicles

This is the elution section of extracellular vesicles from **Washed EVs-binding beads**.

- 1) Add 50 μ L of Exosome Elution Buffer into the 1.5mL Reaction Tube containing **Washed EVs-binding beads** and mix by vortexing^{※13}.
 - 2) After spin-down of the 1.5mL Reaction Tube from step 1), place it on the magnetic stand for 1 minute to collect the beads. Transfer the supernatant into a new 1.5mL microcentrifuge tube (not included in this kit) after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube.
 - 3) Add 50 μ L of Exosome Elution Buffer into the 1.5mL Reaction Tube from step 2) again, and mix by vortexing^{※13}.
 - 4) After spin-down of the 1.5mL Reaction Tube from step 3), place it on the magnetic stand for 1 minute to collect the beads. Transfer the supernatant into the 1.5mL microcentrifuge tube from step 2) and pool together^{※14,15} (Total: 100 μ L) after confirming attachment of the beads to the wall surface of the 1.5mL Reaction Tube.

〔Optional〕


In the case of RNA purification from eluted extracellular vesicles, please add 100 μ L of Exosome Elution Buffer to pooled eluate. After preparing 200 μ L of sample, RNAs are efficiently extracted from extracellular vesicles by using Wako's microRNA Extractor SP kit (Code No. 295-71701) according to the manufacturer's instruction manual.

- ※13 When the aggregates of magnetic beads are remaining, please suspend magnetic beads uniformly by tapping and/or pipetting.

- ※14 **Exosome Capture-immobilized beads** are reusable after eluting EVs (up to 5 times). Therefore, when you need to recover more remaining EVs from samples, please repeat the section 3-5 to increase the recovered amount. After recovering the eluate, please add the stored supernatant of section 3 (refer to ※6) into the 1.5mL Reaction Tube containing **Exosome Capture-immobilized beads** from step 4) of section 5 and repeat the section 3-5 as necessary. This kit includes the required reagents to reuse a maximum of 5 times.

- ※15 There is a possibility that pooled eluate is contaminated with trace amount of magnetic beads. When extracellular vesicles are analyzed by Nanoparticle Tracking Analysis (NTA) or electron microscopic analysis, please filter pooled eluate by using 0.45 μ m pore size filter unit (Example: Ultrafree[®] MC 0.45 μ m, Code No. UFC30HV25, Millipore[®]) before analysis as needed.

[Flowchart of Affinity Purification (Section 3-5)]

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Code No.	Description	Size
295-71701	microRNA Extractor SP Kit	50 reactions

Wako Pure Chemical Industries, Ltd.

1-2, Doshimachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan

Telephone : +81-6-6203-3741

Facsimile : +81-6-6201-5964

<http://www.wako-chem.co.jp>

Wako Chemicals USA, Inc.

1600 Bellwood Road
Richmond, VA 23237
U.S.A.

Telephone : +1-804-271-7677

Facsimile : +1-804-271-7791

<http://www.wakousa.com>

Wako Chemicals GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : +49-2131-311-0
Facsimile : +49-2131-311100
<http://www.wako-chemicals.de>

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MagCapture™ Exosome Isolation Kit PS

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