



Design for you

MGISP-100 Compatible Kits

Automated Sample Preparation System





- Easy to use, automated 90% steps, eliminate operators from repetitive processing procedures.
- Compact design, integrated liquid handling, automat– ed magnetic frame and PCR machine into less than 0.5m³.






Safe & Stable

- Integrated UV lamp and HEPA system, verified most MGI library preparation kit.
- Support customized development accord ing to different needs and applications



Product Function



Compatible Kits List of MGI

Kit Name	Item No.
MGIEasy FS DNA Library Prep Set	1000006987
MGIEasy Universal DNA Library Prep Set	1000006985
MGIEasy FS PCR-Free Library Prep Set	1000013454
MGIEasy PCR-Free Library Prep Set	1000013452
MGIEasy Exome FS Library Prep set	1000009658
MGIEasy Cell-free DNA Library Prep Set	1000007037
MGIEasy RNA Library Rapid Prep Set	1000020509
MGIEasy RNA Library Prep Set	1000006383
MGIEasy rRNA Depletion Kit	1000005953
MGIEasy Small RNA Library Prep Kit	1000005269
MGIEasy RNA Directional Library Prep Set	1000006385
MGICare BRCA1/2 Sequencing Library Preparation Kit	1000009522
MGICare Chromosome Copy Number Variation Sequencing Library Prep Set	1000005290
MGICare Detection Kit for Single Cell Chromosome Copy Number Variation Test	1000005291
MGIEasy Universal Library Conversion Kit (App-A)	1000004155
MGIEasy Circularization Kit	1000005259
MGIEasy Magnetic Beads Genomic DNA Extraction Kit	1000010524
MGIEasy Nucleic Acid Extraction Kit (MGISP-100B)	1000020471

Compatible Kits List of Third-party

Kit Brand	Kit Name	Item No.
Agilent Tech	Agilent SureSelect XT V6	5190-8863
10x Genomics	Chromium Next GEM Single Cell 3' Reagent Kits v3.1	1000268
	Chromium™ Single Cell V(D)J Enrichment Kit v3.1	1000005
Pillar	ONCO/Reveal PillarHS Lung Panel	HLA-LL-1001-24
Vazyme	VAHTS Universal Plus DNA Library Prep Kit for Illumina	ND-617
	VAHTS Universal Plus DNA Library Prep Kit for MGI	NDM617
Yeason	Hieff NGS Ultima DNA Library Prep Kit for Illumina	12199ES24/12199ES96
	Hieff NGS OnePot □ DNA Library Prep Kit for MGI	13321ES16
SOPHiA	SOPHIA SOLID TUMOR SOLUTION	\

Throughput of MGISP-100 is 16 samples/run

