

Terasaki HLA Tissue Typing Tray

Key Benefits

- Broad tray selection
- 72 well format
- Accurate typing assignments
- Well defined antisera
- Long-lasting lots assure consistent results
- Superior technical support
- Prescreened rabbit complement provided at no charge



Easily, accurately and quickly assign HLA Class I and Class II Antigens

Terasaki HLA Tissue Typing Trays are used to assign HLA Class I or Class II antigens accurately and quickly. The trays contain positive and negative controls and are designed for use with both staining techniques - dye exclusion and fluorescence. Each well of the tray contains 1 μ l of specific antiserum and 5 μ l of heavy mineral oil. All trays come with pre-screened rabbit complement.

Tissue Typing Trays can be used for a variety of HLA purposes including:

- Solid organ transplantation
- Bone marrow screening
- Disease association
- Cell panel typings
- Paternity testing
- Platelet matching

Tissue Typing Trays


HLA Class I Serology Antigens Detected

A1	B7	B57 (B17)	CW1
A2	B705*	B58 (B17)	CW2
A3	B8	B59	CW4
A11	B13	B60 (B40)	CW5
A1101*	B18	B61 (B40)	CW6
A1102*	B27	B62 (B15)	CW7
A23 (A9)	B2708*	B63 (B15)	CW8
A24 (A9)	B35	B64 (B14)	CW9 (CW3)
A2403 (A9)	B37	B65 (B14)	CW10 (CW3)
A25 (A10)	B38 (B16)	B67	CW15
A26 (A10)	B39 (B16)	B71 (B15)	CW17
A29 (A19)	B4005 (B40)	B72 (B15)	
A30 (A19)	B41	B73	
A31 (A19)	B42	B75 (B15)	
A32 (A19)	B44 (B12)	B76 (B15)	
A33 (A19)	B45 (B12)	B77 (B15)	
A3401 (A10)	B46	B78	
A3401 (A10)	B47	B81	
A36	B48	B8201	
A6601 (A10)	B49 (B21)		
A6602 (A10)	B50 (B21)		
A68 (A28)	B51 (B5)		
A69 (A28)	B5102 (B5)		
A74	B52 (B5)		
A80	B53		
	B54 (B22)		
	B55 (B22)		
	B56 (B22)		

*Supplemental typing trays

LMT144A/B Lot#4, AS72 Lot#006, BS72 Lot#006, BS740-72 Lot#005, C722 Lot#012

Ordering Information

Product	Cat. No
For In Vitro Diagnostic Use. (Unless otherwise stated.) 	
Terasaki HLA-ABC First Tissue Typing Tray	T1-72
Terasaki HLA-ABC Second Tissue Typing Tray	T2-72
Terasaki HLA-ABC Third Tissue Typing Tray	T3-72
Terasaki HLA-C Locus Tissue Typing Tray	C72
Terasaki HLA-C Locus Supplement Tissue Typing Tray	C722
Terasaki HLA-ABC Oriental Tissue Typing Tray	OT-72
Terasaki HLA-A Locus Supplement Tissue Typing Tray	AS72
Terasaki HLA-B Locus Supplement Tissue Typing Tray	BS72
Terasaki HLA-B7/40 Supplement Tissue Typing Tray	BS740-72
Terasaki HLA Class I Tissue Typing Tray Set	TM144
Terasaki HLA Class II Tissue Typing Tray	DR72

Find out more at onelambda.com

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

21001 Kittridge Street
 Canoga Park, CA 91303-2801, USA
 TEL: 818.702.0042 | 800.822.8824
 (except greater LA area)
 FAX: 818.702.6904 | 800.992.2111
 (U.S. and Canada Only)
 INTERNATIONAL: Contact your local distributor



HLA Class II Serology Antigens Detected

DR1	Dr51
DR15 (DR2)	Dr52
DR16 (DR2)	Dr53
DR17 (DR3)	DQ5 (DQ1)
DR18 (DR3)	DQ6 (DQ1)
DR4	DQ2
DR7	DQ7 (DQ3)
DR8	DQ8/9 (DQ3)
DR9	DQ4
DR10	
DR11	
DR12	
DR13 (Rr6)	
DR14 (DR6)	
DR1404(DR6)	
A33 (A19)	
DR103	

(MRD172 LOT#10, LOT#11)

 **ONE LAMBDA**
 A Thermo Fisher Scientific Brand

MKT-PS-TTT-000 | MAY 2016 | REV 0